



**AirFit™**  
**F40**

Actual Size: 148mm (Width)  
105mm (Height)

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size




**Rx** Only



**ResMed Pty Ltd**  
1 Elizabeth Macarthur Drive  
Bella Vista NSW 2153  
Australia

AirFit is a trademark and/or registered trademark of the ResMed family of companies. For patent and other intellectual property information, see [ResMed.com/ip](https://www.resmed.com/ip). For symbols glossary see [ResMed.com/symbols](https://www.resmed.com/symbols). See [ResMed.com](https://www.resmed.com) for other ResMed locations worldwide.  
© 2024 ResMed. 1036419/1 2024-02



Sample size check

10 mm 20 mm 30 mm 40 mm 50 mm 60 mm 70 mm 80 mm 90 mm 100 mm 110 mm 120 mm 130 mm 140 mm 150 mm 160 mm 170 mm

Ruler is for scaling purposes.